

Multi-Purpose Solution



WE CARE MULTI-PURPOSE SOLUTION

Intended Use:

Intended for use with WE CARE Multi-Purpose Solution

Size	GTIN	REF
3.4 fl oz	4049584037794	10048393
12 fl oz	404958403800	10048394

PLEASE READ CAREFULLY AND KEEP THIS PACKAGE INSERT FOR FURTHER USE:

WE CARE Multi-Purpose Contact Lens Solution (WE CARE MPS) cleans, removes protein, disinfects, rinses, and stores soft (hydrophilic) contact lenses including silicone hydrogel (SH) contact lenses as well as rigid gas permeable (RGP) contact lenses including fluoro-silicone acrylate (FSA) and silicone acrylate (SA) contact lenses. This solution **conditions** only RGP contact lenses, including FSA and SA contact lenses.

CONTENTS: A sterile isotonic solution that contains poloxamer, sodium phosphate buffer, sodium chloride, and disodium edetate, preserved with polyhexanide 0.0001%. Contains no chlorhexidine or thimerosal.

ACTIONS: Cleans, loosens and removes accumulations of films, deposits and debris from soft and RGP contact lenses, including SH, FSA and SA contact lenses. Destroys harmful microorganisms on the surface of the lens. Disinfects, rinses, stores lenses before lenses are placed on the eyes. This solution conditions only RGP contact lenses including FSA and SA contact lenses.

INDICATIONS (Uses): **WE CARE MPS** is indicated for use in chemical (not heat) disinfection, daily cleaning, rinsing and storage of daily and extended wear soft contact lenses including SH, and RGP contact lenses including FSA and SA contact lenses as recommended by your eye care practitioner.

CONTRAINDICATIONS (Reasons not to use): If you are allergic to any ingredient in the **WE CARE MPS**, DO NOT use this product.

WARNINGS: To avoid contamination, never touch dropper tip of the container with your hands or to any surface. Keep bottle tightly closed when not in use. Contact lenses and solutions, if not handled properly, can become contaminated by common bacteria known to cause ocular infections. Serious damage to the eye and loss of vision may result from using contaminated solutions or from wearing contaminated or damaged contact lenses.

PROBLEMS WITH CONTACT LENSES AND LENS CARE PRODUCTS COULD RESULT IN SERIOUS INJURY TO THE EYE.

It is essential that you follow your eye care practitioner's directions and all labeling instructions for proper use and care of your lenses and lens care products, including the lens case. EYE PROBLEMS, INCLUDING CORNEAL ULCERS, CAN DEVELOP RAPIDLY AND LEAD TO LOSS OF VISION. IF YOU EXPERIENCE EYE DISCOMFORT, EXCESSIVE TEARING, VISION CHANGES, REDNESS OF THE EYE, IMMEDIATELY REMOVE YOUR LENSES AND PROMPTLY CONTACT YOUR EYE CARE PRACTITIONER.

Daily wear lenses are not indicated for overnight wear and should not be worn while sleeping. Clinical studies have shown the risk of serious adverse reactions is increased when these lenses are worn overnight. Extended wear lenses should be regularly removed for cleaning and disinfection or for disposal and replacement on the schedule prescribed by your eye care practitioner. Clinical studies have shown that there is an increased incidence of serious adverse reactions in extended wear contact lens users as compared to daily wear contact lens users.

Studies have also shown that the risk of serious adverse reactions increases the longer extended wear lenses are worn before removal for cleaning and disinfection or for disposal and replacement.

Clinical studies have also shown that smokers have a higher incidence of adverse reactions.

If you experience eye discomfort, excessive tearing, vision changes, redness of the eye, immediately remove your lenses and promptly contact your eye care practitioner. All contact lens wearers should see their eye care practitioner as directed.

- To avoid contamination, do not touch tip of container to any surface. Replace cap after using.
- To avoid contaminating your solution, do not transfer solution to other bottles or containers.
- Do not use tap water, bottled water, or saliva with lenses or lens case.
- Not for use with heat (thermal) disinfection.
- Store at room temperature.
- Avoid freezing.
- Keep out of reach of children.

PRECAUTIONS:

- Never reuse this solution.
- Always discard the remaining solution from the lens case.
- Always use fresh solution for soaking and storing lenses.
- Always wash your hands before handling your lenses.
- Keep the bottle tightly closed when not in use.
- After inserting your lenses, always thoroughly clean the interior of the lens case with fresh **WE CARE MPS** and air dry after each use.
- Use before the expiration date marked on the carton and bottle.
- Always rinse lenses with **WE CARE MPS** after cleaning to eliminate loosened surface deposits.
- Discard contents of the bottle three (3) months after opening.

ADVERSE REACTIONS (Problems and what to do):

The following problems may occur: eyes sting, burn or itch (irritation), comfort is less than when lens was first placed on the eye, feeling of something in the eye (foreign body, scratched area), excessive watering (tearing) of the eye, unusual eye secretions, redness of the eye, reduced sharpness of vision (poor visual acuity), blurred vision, rainbows or halos around objects, sensitivity to light (photophobia), or dry eyes.

Multi-Purpose Solution



If you notice any of the above:

IMMEDIATELY REMOVE YOUR LENSES:

- If the discomfort or problem stops, then look closely at the lens.
- If the lens is in any way damaged, DO NOT put the lens back on your eye. Place the lens in the storage case and contact your eye care practitioner.
- If the lens has dirt, an eyelash, or other foreign body on it, or the problem stops and the lens appears undamaged, thoroughly clean, rinse and disinfect the lens, then reinsert it.

If the problem continues, IMMEDIATELY remove the lens and consult your eye care practitioner.

If any of the above symptoms occur, a serious condition such as infection, corneal ulcer, neovascularization or iritis may be present. Seek immediate professional identification of the problem and prompt treatment to avoid serious eye damage. All adverse reactions observed while using should be reported to:

WAVE Eye Solutions™, phone: +1 855 728 2889.

DIRECTIONS: To open bottle, turn cap counterclockwise. The cap has a tamper evident ring that breaks when it is turned counterclockwise. DO NOT use if the ring is broken before the first use. Replace cap after use. To obtain cleaning and disinfection these steps must be completed:

STEP 1—Clean Daily: Place 3 drops of **WE CARE MPS** on each side of contact lens. Gently rub both sides of the lens with tip of forefinger for 15-20 seconds. NO SEPARATE DAILY CLEANER IS REQUIRED.

STEP 2—Rinse: Remove surface debris by rinsing thoroughly with **WE CARE MPS**.

STEP 3—Disinfect/Store: Place rinsed contact lenses in lens case and fill with fresh **WE CARE MPS**. Soak for at least six (6) hours (or overnight) before wearing. Only when you follow all the above steps are your lenses ready to wear. No saline rinse is necessary. If debris remains on contact lens, rinse with **WE CARE MPS**. If not wearing contact lenses immediately, store in unopened lens case. Do not store your lenses in ordinary saline in place of **WE CARE MPS**. Saline solution will not disinfect. Re-disinfect after 7 days or more of storage.

STEP 4—Conditioning: ONLY FOR RGP CONTACT LENSES INCLUDING FSA AND SA CONTACT LENSES:

To improve comfort and wettability a few drops of **WE CARE MPS** may be applied to both surfaces of the lens prior to insertion on to the eye.

GOOD LENS CARE PRACTICES:

- Always wash and rinse your hands before you handle your lenses. This will help to prevent eye infections by removing dirt and oils that could get on the lenses.
- Use only the solutions recommended by your eye care practitioner. Seek advice from your eye care practitioner before making any changes to your care regimen to ensure compatibility with lenses.
- Always follow the directions for use in the labeling included with each solution as instructions may be different for each solution.
- Clean, rinse and disinfect your lenses each time you remove them.
- Always handle the same lens, the left or the right, first in order to avoid mix-ups.
- After each use, always empty and rinse the lens case with fresh **WE CARE MPS**, and allow to air dry.
- Lens cases should be replaced frequently to avoid potential microbial contamination and ocular irritation.

HOW SUPPLIED:

WE CARE MPS is available for purchase in a 3.4 fl. oz. (100 mL) or 12 fl. oz (360 mL) plastic bottle and 1 lens case. Bottles and cartons are marked on bottom with a lot number and expiration date.

TAMPER RESISTANT:

Product is filled in a bottle with a tamper evident ring that breaks when the cap is turned counterclockwise. DO NOT use if the ring is broken.

Manufactured in the Netherlands

Distributed by: OCULUS Surgical Inc., Port St. Lucie, FL 34987

WAVE Eye Solutions™ is a Trademark of OCULUS GmbH

www.waveeyesolutions.com

Call: +1 855 728 2889 | GTIN: 4049584037794

IFU_WECARE_
PDF_20260820_rev1



STERILE

LOT

See base
for EXP. date.



RECYCLABLE